Stripette[™] Serological Pipets







Paper/plastic wrap



Clear plastic wrap



Bulk pack

- Manufactured and packaged with animal-free materials
- Sterile
- ▶ RNase-/DNase-free
- Nonpyrogenic
- Noncytotoxic
- Exclusive anti-drip tip assures accurate delivery.
- Color-coded magnifier stripes make volume reading and size selection easier.
- Bidirectional graduations provide choice of ascending and descending scales.
- Negative graduations allow additional working volume.
- ▶ Three packaging options:
 - Individually wrapped, paper/plastic
 - Individually wrapped, clear plastic
 - Bulk packed
- ▶ Volumetric accuracy of ±2% at stated capacity

Individually Wrapped, Paper/Plastic Wrap

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color-coded Stripe	Qty/Pk	Qty/Cs
10074450	1	1/100	0.2	Yellow	50	1,000
10084450	2	1/100	0.2	Green	50	1,000
10420201	5	1/10	2.5	Blue	50	200
10084450	10	1/10	3.0	Orange	50	200
10420201	25	2/10	10.0	Red	25	200
10041591	50	1/2	10.0	Purple	25	100
10420201	100	1	n/a	Aqua	10	100
Individually	Wrapped, Cl	ear Plastic Wrap				
10238800	1	1/100	0.2	Yellow	100	1,000
10265471	1	1/100	0.2	Yellow	100	200
10340501	2	1/100	0.2	Green	100	1,000
10127400	5	1/10	2.5	Blue	50	200
10677341	10	1/10	3.0	Orange	50	200
10206391*	10	1/10	3.0	Orange	50	200
10732742	25	2/10	10.0	Red	50	200
10636391	50	1/2	10.0	Purple	25	100
10194331	100	1	n/a	Aqua	10	100
Bulk Packed	l in Bags					
10420672	1	1/100	0.2	Yellow	50	1,000
10574712	2	1/100	0.2	Green	50	1,000
10758782	5	1/10	2.5	Blue	50	500
10275471	10	1/10	3.0	Orange	50	500
10633053	25	2/10	10.0	Red	25	200
10633053	50	1/2	10.0	Purple	25	100

*Cat. No. 10206391 features a wide tip for handling viscous fluids.





Triple-bagged packaging

Triple-bagged Packaging for Stripette™ Serological Pipets

Corning's triple-bagged packaged pipets are a safety and convenience revolution. This extraordinary packaging is ideal for cell-based vaccine and drug production facilities, or other laboratory clean rooms.

- ▶ Triple-bagged
- Each pipet is individually wrapped in easy open paper/plastic
- No need for decontamination spraying with alcohol
- Sterility Assurance Level (SAL) of 10⁻⁶

Note: Triple-bagged packaged Stripettes have the same features as the standard Stripettes except for the sterility and packaging.









Contained in a twist-tied liner (Bag 1), the entire case contents of triple-bagged pipets can easily be removed from the corrugated box for clean area storage.

At the first stage clean room, the liner (Bag 1) is opened and the double-bagged inner packs with their outer sterile barriers (Bag 2) are removed.

At the second stage clean room the outer sterile barrier (Bag 2) is opened and the inner sterile barrier (Bag 3) is removed.

Reaching the point of use, the inner sterile barrier (Bag 3) is opened and the individually wrapped pipets are removed.

Ordering Information

Triple-bagged Packaging, Individually Wrapped, Paper/Plastic

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color-coded Stripe	Qty/Pk	Qty/Cs
10221303	1	1/100	0.2	Yellow	50	1,000
10231303	2	1/100	0.2	Green	50	1,000
11333715	5	1/10	2.5	Blue	50	200
11333715	10	1/10	3.0	Orange	50	200
10445243	25	2/10	10.0	Red	25	200
10611604	50	1/2	10.0	Purple	25	100
10417163	100	1	n/a	Aqua	10	100

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

FALCON

AXYG≣N

GOSSELIN

PYREX

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

Sweden: fishersci.se Switzerland: fishersci.ch UK: fishers
© 2022 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.

